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RÉFÉRENCES

- (1) Delsys Inc (2014). *Trigno Wireless System: User's Guide*. Natick, MA: Author.
- (2) Delsys Inc (2016). *EMGworks 4 User's Guide*. Retrieved from
<http://www.delsys.com/KnowledgeCenter/NetHelp/>.
- (3) Christiano, M. (2015, July 31). Transformer Isolation. Retrieved December 31, 2017, from
<https://www.allaboutcircuits.com/technical-articles/transformer-isolation/>
- (4) Criswell, E., & Cram, J. R. (2011). *Cram's introduction to surface electromyography*. Sudbury, MA: Jones & Bartlett.
- (5) Delsys Inc. (n.d.). *Technical Note 101: EMG Sensor Placement (Tech.)*. Retrieved October 15, 2017, from <https://www.delsys.com/.../TN101%20-%20EMG%20Sensor%20Placement-web.pdf>
- (6) Raison, M. (n.d.). *Évaluation d'un nouvel instrument de mesure totalement non invasif pour caractériser le comportement musculaire d'enfants atteints d'un trouble de motricité de type DMC ou MNM dans le cadre de leur réadaptation* (Protocole de recherche).
- (7) Sarcher, A. (2014). *Mouvements du membre supérieur chez les enfants hémiplégiques: analyse cinématique et électromyographie* (masters dissertation). Université de Montréal.
- (8) Sarcher, A., Raison, M., Ballaz, L., Lemay, M., Leboeuf, F., Trudel, K., & Mathieu, P. (2015). *Impact of muscle activation on ranges of motion during active elbow movement in children with spastic hemiplegic cerebral palsy*. Clinical Biomechanics, 30(1), 86-94.
doi:10.1016/j.clinbiomech.2014.10.009
- (9) Rose, W. (n.d.). *Electromyogram Analysis (Tech.)*. Retrieved July 2017, 10, from
<https://www1.udel.edu/biology/rosewc/kaap686/notes/EMG%20analysis.docx>
- (10) Anatomy: Muscle Atlas. (n.d.). Retrieved November 29, 2017, from
<https://www.kenhub.com/en/dashboard>
- (11) Brachioradialis. (n.d.). Retrieved November 29, 2017, from
<http://www.exrx.net/Muscles/Brachioradialis.html>
- (12) Gaudet G, Raison M, Achiche S. (2018) Classification of Upper limb phantom movements in transhumeral amputees using electromyographic and kinematic features. Engineering Applications of Artificial Intelligence 68, 153-164. <https://doi.org/10.1016/j.engappai.2017.10.017>
- (13) J Clément, M Raison, DM Rouleau. (2017) Reproducibility analysis of upper limbs reachable workspace, and effects of acquisition protocol, sex and hand dominancy. Journal of biomechanics. <https://doi.org/10.1016/j.jbiomech.2017.12.010>
- (14) B Samadi, M Raison, L Ballaz, S Achiche. (2017) Decomposition of three-dimensional ground-reaction forces under both feet during gait. Journal of musculoskeletal & neuronal interactions 17(4), 283-291

- (15) Dal Maso F, Begon M, Raison M. (2017) Methodology to Customize Maximal Isometric Forces for Hill-Type Muscle Models. *Journal of applied biomechanics* 33 (1), 80-86.
<https://doi.org/10.1123/jab.2016-0062>
- (16) M Leroux, S Achiche, M Raison. (2017). Assessment of accuracy for target detection in 3D-space using eye tracking and computer vision. *PeerJ Preprints* 5, e2718v1.
<https://doi.org/10.7287/peerj.preprints.2718v1>
- (17) Raouafi, S., Achiche, S., Begon, M., Sarcher, A. & Raison, M. (2017). Classification of upper limb disability levels of children with spastic unilateral cerebral palsy using K-means algorithm. *Medical and Biological Engineering and Computing*, 11 pages. Tiré de <https://doi.org/10.1007/s11517-017-1678-y>
- (18) Zizoua, C., Raison, M., Boukhenous, S., Attari, M. & Achiche, S. (2017). Development of a Bracelet with Strain-Gauge Matrix for Movement Intention Identification in Traumatic Amputees. *IEEE Sensors Journal*, 17(8), 2464-2471. Tiré de <https://doi.org/10.1109/JSEN.2017.2666784>
- (19) Bousquet-Jette, C., Achiche, S., Beaini, D., Cio, Y., Leblond-Menard, C. & Raison, M. (2017). Fast scene analysis using vision and artificial intelligence for object prehension by an assistive robot. *Engineering Applications of Artificial Intelligence*, 63, 33-44. Tiré de
<https://doi.org/10.1016/j.engappai.2017.04.015>
- (20) Sarcher, A., Raison, M., Leboeuf, F., Perrouin-Verbe, B., Brochard, S. & Gross, R. (2017). Pathological and physiological muscle co-activation during active elbow extension in children with unilateral cerebral palsy. *Clinical Neurophysiology*, 128(1), 4-13. Tiré de <https://doi.org/10.1016/j.clinph.2016.10.086>
- (21) Clement, J., Menard, J., Raison, M., Dumais, J., Dubois, L. & Rouleau, D.M. (2017). Three-dimensional analysis of the locked position in patients with recurrent shoulder instability. *Journal of Shoulder Elbow Surgery*, 26(3), 536-543. Tiré de <https://doi.org/10.1016/j.jse.2016.07.031>
- (22) Achiche, S., Shlechtingen, M., Raison, M., Baron, L. & Santos, I.F. (2016). Adaptive Neuro-Fuzzy Inference System Models for Force Prediction of a Mechatronic Flexible Structure. *Journal of Integrated Design & Process Science*, 19(3), 77-94. Tiré de <https://doi.org/10.3233/jid-2015-0017>
- (23) Veilleux, L.N., Raison, M., Rauch, F., Robert, M. & Ballaz, L. (2016). Agreement of spatio-temporal gait parameters between a vertical ground reaction force decomposition algorithm and a motion capture system. *Gait & Posture*, 43, 257-264. Tiré de <https://doi.org/10.1016/j.gaitpost.2015.10.007>
- (24) Samadi, B., Achiche, S., Parent, A., Ballaz, L., Chouinard, U. & Raison, M. (2016). Custom sizing of lower limb exoskeleton actuators using gait dynamic modelling of children with cerebral palsy. *Computer Methods in Biomechanics and Biomedical Engineering*, 19(14), 1519-1524. Tiré de
<https://doi.org/10.1080/10255842.2016.1159678>
- (25) Zizoua, C., Raison, M., Boukhenous, S., Attari, M. & Achiche, S. (2016). Detecting muscle contractions using strain gauges. *Electronics Letters*, 52(22), 1836-1838. Tiré de
<https://doi.org/10.1049/el.2016.2986>

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- (37) Dal Maso, F., Blache, Y., Raison, M., Lundberg, A. & Begon, M. (2015). Glenohumeral joint kinematics measured by intracortical pins, reflective markers, and computed tomography: A novel technique to assess acromiohumeral distance. *Journal of Electromyography and Kinesiology*, 29, 4-11. Tiré de <https://doi.org/10.1016/j.jelekin.2015.07.008>
- (38) Dal Maso, F., Raison, M., Lundberg, A., Arndt, A., Allard, P. & Begon, M. (2015). Glenohumeral translations during range-of-motion movements, activities of daily living, and sports activities in healthy participants. *Clinical Biomechanics*, 30(9), 1002-1007. Tiré de <https://doi.org/10.1016/j.clinbiomech.2015.06.016>
- (39) Sarcher, A., Ballaz, L., Leboeuf, F., Gross, R., Perrouin-Verbe, B. & Raison, M. (2015). Impact of elbow angular velocity on muscle activation and coactivation during active elbow extension and supination in children with spastic hemiplegic cerebral palsy. *Gait & Posture*, 42(S1), S18-S19. Tiré de <https://doi.org/10.1016/j.gaitpost.2015.06.039>
- (40) Haering, D., Blache, Y., Raison, M. & Begon, M. (2015). Mechanical risk of rotator cuff repair failure during passive movements: A simulation-based study. *Clinical Biomechanics*, 30(10), 1181-1188. Tiré de <https://doi.org/10.1016/j.clinbiomech.2015.08.006>
- (41) Laitenberger, M., Raison, M., Perie, D. & Begon, M. (2015). Refinement of the upper limb joint kinematics and dynamics using a subject-specific closed-loop forearm model. *Multibody System Dynamics*, 33(4), 413-438. Tiré de <https://doi.org/10.1007/s11044-014-9421-z>
- (42) Dal Maso, F., Raison, M., Lundberg, A., Arndt, A. & Begon, M. (2014). Coupling between 3D displacements and rotations at the glenohumeral joint during dynamic tasks in healthy participants. *Clinical Biomechanics*, 29(9), 1048-1055. Tiré de <https://doi.org/10.1016/j.clinbiomech.2014.08.006>
- (43) Mahaudens, P., Raison, M., Banse, X., Mousny, M. & Detrembleur, C. (2014). Effect of long-term orthotic treatment on gait biomechanics in adolescent idiopathic scoliosis. *Spine Journal*, 14(8), 1510-1519. Tiré de <https://doi.org/10.1016/j.spinee.2013.08.050>
- (44) Hart, R., Ballaz, L., Robert, M., Pouliot, A., D'Arcy, S., Raison, M. & Lemay, M. (2014). Impact of Exercise-Induced Fatigue on the Strength, Postural Control, and Gait of Children with a Neuromuscular Disease. *American Journal of Physical Medicine & Rehabilitation*, 93(8), 649-655. Tiré de <https://doi.org/10.1097/phm.0000000000000091>
- (45) Mukherjee, A., Raison, M., Sahni, T., Arya, A., Lambert, J., Marois, P., James, P., Parent, A. & Ballaz, L. (2014). Intensive rehabilitation combined with HBO therapy in children with cerebral palsy: a controlled longitudinal study. *Undersea and Hyperbaric Medicine Journal*, 41(2), 77-85.
- (46) Haering, D., Raison, M., Arndt, A. & Begon, M. (2014). Kinematic model and elbow flexion interaction on shoulder range of motion. *Computer Methods in Biomechanics and Biomedical Engineering*, 17(SUPP1), 84-85. Tiré de <https://doi.org/10.1080/10255842.2014.931150>
- (47) Haering, D., Raison, M. & Begon, M. (2014). Measurement and description of three-dimensional shoulder range of motion with degrees of freedom interactions. *Journal of Biomechanical Engineering*, 136(8). Tiré de <https://doi.org/10.1115/1.4027665>

- (48) Ballaz, L., Raison, M. & Detrembleur, C. (2013). Decomposition of the Vertical Ground Reaction Forces During Gait on a Single Force Plate. *Journal of Musculoskeletal & Neuronal Interactions*, 13(2), 236-243.
- (49) Mahaudens, P., Raison, M., Banse, X., Mousny, M. & Detrembleur, C. (2013). Effects of long-term orthotic treatment on gait biomechanics in adolescent idiopathic scoliosis. *The Spine Journal*.
- (50) Letellier, G., Parent, A., Raison, M., Marois, P., Mercier, C., Crevier, L., Laberge, M. & Ballaz, L. (2013). Twenty-five years of selective dorsal rhizotomy based on clinical observations instead of intraoperative electrophysiological monitoring: The Sainte-Justine hospital experience. *Annals of Physical and Rehabilitation Medicine*, 56(S1), e301. Tiré de <https://doi.org/10.1016/j.rehab.2013.07.772>
- (51) Mahaudens, P., Banse, X., Mousny, M., Raison, M. & Detrembleur, C. (2013). Very Short-Term Effect of Brace Wearing on Gait in Adolescent Idiopathic Scoliosis Girls. *European Spine Journal*, 22(11), 2399-2406.
- (52) Ramananarivo, M., Raison, M., Barron, O. & Achiche, S. (2017). Biofidelic design of the forearm of a myoelectric prosthesis with maximum functional volume. Communication présentée à 21st International Conference on Engineering Design (ICED 2017), Vancouver, BC, Canada (p. 479-488).
- (53) Leroux, M., Achiche, S. & Raison, M. (2017). Design guidelines for shoulder design of an anthropomorphic robotic ARM. Communication présentée à 21st International Conference on Engineering Design (ICED 2017), Vancouver, BC, Canada (p. 41-50).
- (54) Zizoua, C., Achiche, S., Boukhenous, S., Attari, M. & Raison, M. (2017). Wearable bio-instrument for detecting body movement based on strain gauges. Communication présentée à 30th IEEE Canadian Conference on Electrical and Computer Engineering (CCECE 2017), Windsor, ON, Canada (4 pages). Tiré de <https://doi.org/10.1109/CCECE.2017.7946640>
- (55) Leroux, M., Raison, M., Adadja, T. & Achiche, S. (2015). Combination of eyetracking and computer vision for robotics control. Communication présentée à IEEE Conference on Technologies for Practical Robot Applications, TePRA, Woburn, MA, United states. Tiré de <https://doi.org/10.1109/TePRA.2015.7219692>
- (56) Hernandez, S., Raison, M., Baron, L. & Achiche, S. (2015). Refinement of exoskeleton design using multibody modeling : An overview. Communication présentée à CCToMM Mechanisms, Machines, and Mechatronics Symposium (M3 2015), Ottawa, Ont., Canada. Tiré de http://www.cctomm.mae.carleton.ca/CCToMM_2015.pdf
- (57) Geoffroy, P., Bordron, O., Mansard, N., Raison, M., Stasse, O. & Bretl, T. (2014). A two-stage suboptimal approximation for variable compliance and torque control. Communication présentée à 13th European Control Conference, ECC 2014, Strasbourg, France (p. 1151-1157). Tiré de <https://doi.org/10.1109/ECC.2014.6862557>
- (58) Hernandez, S., Raison, M., Torres, A., Gaudet, G. & Achiche, S. (2014). From on-body sensors to in-body data for health monitoring and medical robotics : a survey. Communication présentée à Global

Information Infrastructure and Networking Symposium (GIIS 2014), Montreal, Quebec, Canada. Tiré de <https://doi.org/10.1109/GIIS.2014.6934279>

(59) Samadi, B., Raison, M., Detrembleur, C. & Ballaz, L. (2014). Real-time detection of reaction forces during gait on a ground equipped with a large force platform. Communication présentée à 2014 Global Information Infrastructure and Networking Symposium, GIIS 2014, Montreal, QC, Canada. Tiré de <https://doi.org/10.1109/GIIS.2014.6934268>

(60) Schlechtingen, M., Achiche, S., Costa, T.L., Raison, M. & Santos, I. (2014). Using data mining approaches for force prediction of a dynamically loaded flexible structure. Communication présentée à 12th Biennial conference on engineering systems design and analysis (ESDA 2014), Copenhagen, Denmark. Tiré de <https://doi.org/10.1115/ESDA2014-20259>

(61) Gervais B, Vadean A, Brochu M, Raison M. (2017). Development of a Tool for the Sensitivity Analysis of Design Parameters of Femoral Implants in the Human Body During Gait. ECCOMAS Thematic Conference on MULTIBODY DYNAMICS. Prague, June 19 - 22, 2017

(62) Sarcher A, Brochard S, Raison M, Leboeuf F, Hug F, Perrouin-Verbe B, Gross R. (2017) Upper limb co-activation in typically developing children and in children with unilateral cerebral palsy: reliability and sensitivity to the EMG normalization technique. *Gait & Posture* 57, 285-286

(63) Hernandez, S., Leroux, M., Achiche, S. & Raison, M. (2016). Based Learning involving physical prototypes to support the teaching of multibody modeling applied to the human body. Communication présentée à Joint International Conference on Multibody System Dynamics (IMSD 2016), Montreal, Quebec.

(64) Gervais, B., Vadean, A., Brochu, M. & Raison, M. (2016). Evolution of the stress distribution of a femoral orthopedic implant by combining multibody dynamics and finite element modeling. Communication présentée à 4th Joint International Conference on Multibody System Dynamics (IMSD 2016), Montréal, Québec, Canada.

(65) Blanchet, L., Dal Maso, F., Begon, M., Achiche, S. & Raison, M. (2016). Impact of muscle modeling refinement on its kinematics and dynamics. Communication présentée à Joint International Conference on Multibody System Dynamics (IMSD 2016), Montreal, Quebec.

(66) Hernandez, S., Leroux, M., Achiche, S. & Raison, M. (2016). Project Based Learning involving physical prototypes to support the teaching of multibody modeling applied to the human body. Communication présentée à Joint International Conference on Multibody System Dynamics (IMSD 2016), Montreal, Quebec.

(67) Larouche-Guilbert, M., Barron, O., Achiche, S. & Raison, M. (2016). Real-time assessment of muscle forces based on multibody dynamics and electromyography. Communication présentée à Joint International Conference on Multibody System Dynamics (IMSD 2016), Montreal, Quebec.

(68) Leroux, M., Beaini, D., Achiche, S. & Raison, M. (2016). Revisiting the Agile eye using Euler-Lagrange Formalism and Torque Control. Communication présentée à Joint International Conference on Multibody System Dynamics (IMSD 2016), Montreal, Quebec.

- (69) Wen, J., Blanchet, L., Dal Maso, F., Achiche, S. & Raison, M. (2016). Towards the identification of the most realistic strategy of muscle force distribution during movement cycles. Communication présentée à Joint International Conference on Multibody System Dynamics (IMSD 2016), Montreal, Quebec.
- (70) Samadi B, Achiche S and Raison M. (2015). Benchmark of the upper limb in 3D to analyze internal effort quantification and realistic movement reproduction. ECCOMAS Thematic Conference on Multibody Dynamics, Barcelona, June 29- July 2, 2015.
- (71) Sarcher, A., Raison, M., Gross, R., Ballaz, L., Leboeuf, F., Brochard, S. & Perrouin-Verbe, B. (2015). Distinction between pathological and functional co-activation during active elbow extension in children with spastic hemiplegic cerebral palsy. Communication présentée à 24th Annual Meeting of the European Society of Movement Analysis in Adults and Children (ESMAC 2015), Heidelberg, Germany.
- (72) Yazji, M., Raison, M., Mahaudens, P., Detrembleur, C., Mousny, M., Labelle, H. & Aubin, C.-É. (2014). Comparison of lower limb joint efforts during gait between healthy and scoliotic adolescents. Communication présentée à 11th World congress on computational mechanics (WCCM 2014), Barcelona, Spain.
- (73) Laitenberger, M., Raison, M., Gaudet, G., Périé, D., Dal Maso, F. & Begon, M. (2014). Forearm modeling refinement : sensitivity on the upper limb osteo-articular kinematics and dynamics. Communication présentée à 11th World congress on computational mechanics (WCCM 2014), Barcelona, Spain.
- (74) Laitenberger, M., Begon, M., Gaudet, G., Périé, D., Dal Maso, F. & Raison, M. (2014). Impact of the forearm modeling refinement on the upper limb joint kinematics and dynamics. Communication présentée à 3rd International multibody system dynamics conference (IMSD 2014), Busan, Korea.
- (75) Gaudet, G., Begon, M., Achiche, S., Dal Maso, F., Musy, G. & Raison, M. (2014). Muscle overactuation problem : two solutions from an undeterminate to a deterministic problem. Communication présentée à 3rd International multibody system dynamics conference (IMSD 2014), Busan, Korea.
- (76) Gaudet, G., Raison, M., Achiche, S., Dal Maso, F., Musy, G. & Begon, M. (2014). Solutions to the muscle redundancy problem : from an undeterminate to a deterministic problem. Communication présentée à 11th World congress on computational mechanics (WCCM 2014), Barcelona, Spain.
- (77) Geoffroy, P., Mansard, N., Raison, M., Achiche, S., Tassa, Y. & Todorov, E. (2014). From inverse kinematics to optimal control, accepted to the advances in robot kinematics. Dans Advances in Robot Kinematics (p. 409-418). Springer. Tiré de https://doi.org/10.1007/978-3-319-06698-1_42
- (78) Raison, M. (2017). Method and system for decomposing superimposed ground reaction forces into individual force profiles (Brevet n° US 9579044). Tiré de <https://www.google.ca/patents/US9579044>
- (79) Raison, M. (2012). Portable and clinical tool for dynamic quantification of the muscles forces performed by the human body (Brevet n° CA EFS-ID-12883045).

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Conditions d'utilisation: Tous droits réservés / All rights reserved